

BD Vacutainer® Venous Blood Collection

Tube Guide For the full array of BD Vacutainer* Blood Collection Tubes, visit www.bd.com/vacutainer. Many are available in a variety of sizes and draw volumes (for pediatric applications). Refer to our website for full descriptions.

BD Vacutainer® Tubes with BD Hemogard [™] Closure	BD Vacutainer [®] Tubes with Conventional Stopper	Additive	Inversions at Blood Collection*	Laboratory Use	Your Lab's Draw Volume/Remarks
Gold	Red/ Gray	Clot activator and gel for serum separation	5	For serum determinations in chemistry. May be used for routine blood donor screening and diagnostic testing of serum for infectious disease." Tube inversions ensure mixing of dot activator with blood Blood dotting time: 30 minutes.	
Light Green	Green/ Gray	 Lithium heparin and gel for plasma separation 	8	For plasma determinations in chemistry. Tube inversions ensure mixing of anticoagulant (heparin) with blood to prevent clotting.	
Red	Red	 Silicone coated (glass) Clot activator, Silicone coated (plastic) 	0 5	For serum determinations in chemistry. May be used for routine blood donor screening and diagnostic testing of serum for infectious disease." Tube inversions ensure mixing of dot a citvator with blood. Blood clotting time: 60 minutes.	
Orange		 Thrombin-based clot activator with gel for serum separation 	5 to 6	For stat serum determinations in chemistry. Tube inversions ensure mixing of dot activator with blood. Blood dotting time: 5 minutes.	
Orange		Thrombin-based clot activator	8	For stat serum determinations in chemistry. Tube inversions ensure mixing of clot activator with blood. Blood dotting time: 5 minutes.	
Royal Blue		 Clot activator (plastic serum) K₂EDTA (plastic) 	8 8	For trace-element, toxicology, and nutritional-chemistry determinations. Special stopper formulation provides low levels of trace elements (see package insert). Tube inversions ensure mixing of either dot activator or anticoagulant (EDTA) with blood.	
Green	Green	• Sodium heparin • Lithium heparin	8 8	For plasma determinations in chemistry. Tube inversions ensure mixing of anticoagulant (heparin) with blood to prevent dotting.	
Gray	Gray	 Potassium oxalate/ sodium fluoride Sodium fluoride/Na₂ EDTA Sodium fluoride (serum tube) 	8 8 8	For glucose determinations. Oxalate and EDTA anticoagulants will give plasma samples. Sodium fluoride is the antiglycolytic agent. Tube inversions ensure proper mixing of additive with blood.	
Tan		• K ₂ EDTA (plastic)	8	For lead determinations. This tube is certified to contain less than .01 µg/mL(ppm) lead Tube inversions prevent clotting.	
	Yellow	 Sodium polyanethol sulfonate (SPS) Acid citrate doxtrose additives (ACD): Solution A - 22.0 g/L trisodium citrate, 8.0 g/L citric acid, 24.5 g/L doxtrose - 13.2 g/L trisodium citrate, 4.8 g/L citric acid, 14.7 g/L doxtrose 	8 8 8	SPS for blood culture specimen collections in microbiology. ACD for use in blood bank studies, HLA phenotyping, and DNA and paternity testing. Tube inversions ensure mixing of anticoagulant with blood to prevent dotting.	
Lavender	Lavender	 Liquid K₃EDTA (glass) Spray-coated K₂EDTA (plastic) 	8 8	K ₂ EDTA and K ₂ EDTA for whole blood hematology determinations. K ₂ EDTA may be used for routine immunohematology testing, and blood donor screening. Tube inversions ensure mixing of anticoagulant (EDTA) with blood to prevent dotting.	
White		 K₂EDTA and gel for plasma separation 	8	For use in molecular diagnostic test methods (such as, but not limited to, polymerase chain reaction [PCR] and/or branched DNA [DDNA] amplification techniques.) Tube inversions ensure mixing of anticoagulant (EDTA) with blood to prevent dotting.	
Pink	Pink	 Spray-coated K₂EDTA (plastic) 	8	For whole blood hematology determinations. May be used for routine immunohematology testing and blood donor screening. ¹⁰ Designed with special cross-match label for patient information required by the AABB. Tube inversions prevent clotting.	
Light Blue Clear	Light Blue	 Buffered sodium citrate 0.105 M (=3.2%) glass 0.109 M (3.2%) plastic Citrate, theophyline, adenosine, dipyridamole (CTAD) 	3-4 3-4	For coagulation determinations. CTAD for selected platelet function assays and routine coagulation determination. Tube inversions ensure mixing of anticoagulant (citrate) to prevent dotting.	
Clear	New Red/ Light Gray	None (plastic)	0	For use as a discard tube or secondary specimen tube.	
Note: BD Vacutaine		ic and partial draw app	ications can l	be found on our website.	

BD Diagnostics Preanalytical Systems 1 Becton Drive Franklin Lakes, NJ 07417 USA

BD Global Technical Services: 1.800.631.0174 BD Customer Service: 1.888.237.2762 www.bd.com/vacutainer

Invert gently, do not shake
 The performance characteristics of these tubes have not been established for infectious disease testing in general; therefore, users must validate the use of these tubes for their specific assay-instrument/nagent system combinations and specimen storage conditions.
 The performance characteristics of these tubes have not been established for immunohematology testing in general; therefore, users must validate the use of these tubes have not been established for immunohematology testing in general; therefore, users must validate the use of these tubes of their specific assay-instrument/regency system combinations and specimen storage conditions.